



**ProTrainings.eu**<sup>®</sup>



# The Dechoker Choking First-Aid Device

## FDA Class II medical device compliant & works in seconds

The Dechoker uses powerful suction to remove food from the airway. It works in seconds and is easy to use on yourself or someone else. Registered as an FDA Class 2 medical device, the Dechoker is approved to be manufactured sterilized and ISO13485 compliant. Dechoker can be used for up to two years in a home. It is designed for single use in public spaces or places of business.

When you consider the safety of your family, friends, workplace, and community, having the Dechoker on hand can give you peace of mind.

**A life-threatening choking attack can happen to anyone, even you.**

The elderly, small children (infants especially) and people suffering from certain diseases are more likely to experience choking. Here are just some of the statistics:

- A leading cause of death in children- one every 5 days.
- A leading cause of death in persons over the age of 65.
- A choking death occurs every two hours.
- Choking causes over 100,000 visits to the A&E yearly.

The team at Dechoker has designed a revolutionary new medical device that can help increase your chances or a loved one's chances of surviving a choking emergency. It helps to know Abdominal Thrusts, but you might not be big or strong enough to perform it on the person choking. The other alternative, performing a Tracheotomy, is something very few would be able to do. Dechoker is our very sensible answer to this problem. Time is of the essence when trying to help someone choking, and it's very distressing to all involved. Dechoker is our very sensible answer to this problem.

The Dechoker is easy to use and can even be used if you are alone. We believe dechoker can give you the piece of mind you need to be there for your family, co-workers and wherever your help might be needed. Every home, school, building and family vehicle should have one!



Available in  
three sizes,  
Adult, child  
and infant

## DECHOKER INSTRUCTIONS

1. Remove Dechoker from package, pull handle once or twice
2. Lay individual on their back, tilt up head, lifting chin up for access to airway
3. Insert tube into mouth, respirator face cap covering mouth and nose for no longer than 3 second intervals
4. Apply thumb at bottom of chin and index finger on one side of respirator, middle finger on the other side of respirator
5. Apply light pressure, begin pulling the plunger upwards. Repeat steps 4-5 if necessary
6. Never leave respirator covering mouth and nose over 3 seconds at any point. Countdown 3...2...1
7. Roll individual over on their side to allow debris out of mouth to avoid pulmonary aspiration

The 'golden hour' is firmly established as a key principle of timely management following serious injury. Seconds really do count when the airway is obstructed and avoidable delay can result in death. According to information provided by the NHS and UK and European Resus Councils, one should NEVER use abdominal thrusts with babies under one year old, pregnant women or people who are obese. Never attempt to put your fingers into your child's throat to feel for the object; doing so could lodge the item more firmly in his airway.

## Time is of The Essence

- 0-4 minutes: Brain Damage Unlikely
- 4-6 minutes: Brain Damage Possible
- 6-10 minutes: Brain Damage Likely
- 10+ minutes: Probable Brain Death
- Average Emergency Response Time: 7-12 minutes

## EXPERT TESTIMONIAL:

*'In my opinion the Dechoker ACD (Anti-Choking Device) should be part of all resuscitation/first-aid kits; in ambulances, for the fire brigade, the police...and in every home. Where abdominal thrusts and back slaps have not worked or the first responder is unable to perform this intervention, the Dechoker ACD should be used ASAP'*

Dr A Biswas, FRCA, Consultant Anaesthetist, Clinical Tutor, Nottingham University NHS Trust, United Kingdom.

## ProTrainings Europe Limited

First Floor, Digby House  
Riverside Office Centre  
Causton Road, Colchester  
Essex CO1 1RJ

Telephone: 01206 805359

Email: [dechoker@protrainings.eu](mailto:dechoker@protrainings.eu)

Web: [www.protrainings.eu](http://www.protrainings.eu)

