## Dog CPR Flowchart



- Open airway
- Check for normal breathing for no more than 10 seconds
- Call for help and get someone to call your vet if possible
- Open airway, pull the tongue to one side and seal the mouth. Deliver 5 rescue breaths until you see the chest rise by blowing in the nose
- Check pulse and if no pulse, deliver 30 compressions in centre of the chest at a rate ot 100-120 compressions per minute to a depth of at least one third the depth of the chest, with one hand
- Then open the airway and deliver 2 breaths until you see the chest rise
- Repeat at 30 compressions and 2 breaths. Check pulse every minute and strop compressions if you detect a pulse
- Continue until they start to breathe, if you are too tired to continue, let someone else take over or the vet arrives



ProTrainings Europe Limited

650 The Crescent • Colchester Business Park, Colchester • Essex CO4 9YQ Telephone: 01206 805359 • Email: info@protrainings.uk

## www.protrainings.uk