
Dog CPR Flowchart



- Open airway
- Check for normal breathing for no more than 10 seconds
- Call for help and get someone to call your vet if possible
- Open airway, pull the tongue to one side and seal the mouth. Deliver 5 rescue breaths until you see the chest rise by blowing in the nose
- Check pulse and if no pulse, deliver 30 compressions in centre of the chest at a rate ot 100-120 compressions per minute to a depth of at least one third the depth of the chest, with one hand
- Then open the airway and deliver 2 breaths until you see the chest rise
- Repeat at 30 compressions and 2 breaths. Check pulse every minute and strop compressions if you detect a pulse
- Continue until they start to breathe, if you are too tired to continue, let someone else take over or the vet arrives



22 Westside Centre, London Road, Stanway, Colchester CO3 8PH

Telephone: 01206 805359 • Email: info@protrainings.uk

